MANUALE D'USO
USER MANUAL
MANUEL D'UTILISATION
BEDIENUNGSANLEITUNG
GEBRUIKERSHANDLEIDING

BKRTR

BKRTS

BKRTA

ARIIN®

www.ariawheels.com

ITA	<i>Manuale d'uso</i>	<i>pag. 5</i>
ENG	<i>User manual</i>	<i>pag.15</i>
FRA	<i>Manuel d'utilisation</i>	<i>pag. 25</i>
DE	<i>Bedienungsanleitung</i>	<i>pag. 35</i>
NL	<i>Gebruikershandleiding</i>	<i>pag. 45</i>

BKRTR

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Legale Rappresentante
(Authorized/Responsible Person)
Firma: Marco Paolucci
Signed

Data: 12/06/2017
Date: 12/06/2017

This user's manual is an integral part of the product and must be supplied with the wheelchair.

ENG

INDEX

1. INTRODUCTION	pg.16
1.1 Information related to the present user's manual	
1.2 Explanation of symbols	
1.3 After sale warranty	
1.4 Overview	
2. SAFETY	pg.17
2.1 Customer transport on a wheelchair equipped with an aria® backrest in a motor vehicle	
2.2 Intended use	
2.3 Accessories warning	
2.4 Installation warnings	
3. CONFIGURATION	pg.19
3.1 Hardware kit contents	
3.2 Tools required	
3.3 Mounting configuration	
4. USAGE	pg.20
4.1 Lumbar support installation/adjustment	
4.2 Lumbar support adjustment	
5. CARE AND MAINTENANCE	pg.21
5.1 Cleaning	
5.2 Inspection	
6. AFTER USE	pg.21
6.1 Reuse	
6.2 Disposal	
7. TECHNICAL DATA	pg.22
7.1 Weight and dimensions	
7.2 Environmental conditions	
7.3 Labels	
8. DECLARATION OF CONFORMITY	pg.24

1 INTRODUCTION

1.1 INFORMATION RELATED TO THE PRESENT USER'S MANUAL

We thank you for choosing the ARIA postural backrest.

This user's manual contains important indications for the use of the product. To be able to utilise the postural backrest correctly, we ask you to read the manual carefully and to observe the safety indications.

In the event of vision problems, it is possible to view this user's manual as a PDF file on the Internet at www.ariawheels.com and zoom in on the screen according to your needs. If you are unable to enlarge the texts and graphics enough, please contact an ARIA official distributor in your national territory; the addresses are listed on the website. If necessary, we will provide a high resolution PDF file of the manual. Additionally, it is possible to obtain the reading of the PDF file using appropriate software that use special linguistic functions installed on your computer.

1.2 EXPLANATION OF SYMBOLS

Symbols in this manual.

All the instructions in this manual are labelled with symbols. In front of the individual messages are symbols and words that indicate the gravity of the danger.

WARNING !

Indicates a potential danger that could cause severe injuries or death, if not avoided.

ATTENTION !

CAUTION !

Indicates a potential danger that could cause small to medium injuries, if not avoided.

WARNING !

IMPORTANT !

Indicates a potential danger that could cause material damages, if not avoided. Indicates tips and useful recommendations for an efficient and inconvenience-free use of the backrest.

IMPORTANT !

This product is in compliance with the directive (93/42/ECC) regarding medical products. The launch date of this product is indicated in the declaration of compliance CE.

SYMBOLS ON THE PRODUCT

An identification label is applied on the product.

The symbols on the label are defined in detail in section 11.4 of this user's manual.

1.3 AFTER SALE WARRANTY

Aria Wheels guarantees the absence of defects and the functionality of its postural backrests. The warranty covers defects that arise from manufacturing, the use of bad quality materials or the inappropriate working of these. The warranty rights towards Aria Wheels S.r.l.s. can only be executed by the seller, and not the user of the product.

The commercial warranty does not cover the normal wear, damages or faults that come directly or indirectly from accidents, falls, collisions, improper use, insufficient maintenance and alterations made by non authorised personnel.

The parts subject to wear and particular consumption are not covered by warranty. The warranty falls if any unexpected alterations are made to the wheelchair, or if alterations are made with inadequate or unoriginal parts. In case of not authorized alterations, the warranty does not cover the costs arising from the elimination of the defect such as shipping costs, labour costs, expenses etc.

The warranty has a duration of 24 (twenty-four) months starting from the date of invoice.

All the components not manufactured by Aria Wheels S.r.l.s. have a warranty from their own manufacturer.

Moreover, the terms and conditions are an integral part of the general and specific conditions for the individual countries in which the product is commercialised.

1.4 OVERVIEW

IMPORTANT INFORMATION

The best way to avoid problems related to pressure sores is to understand their causes and your role in a skin management program. Your therapist and physician should be consulted if you have questions regarding individual limitations and needs.

All back cushions should be selected carefully. Working with your therapist and physician is the best way to assure that a cushion choice matches your individual needs. As the needs of the individual become more complex, the evaluation becomes more important.

ARIA® BACKREST

The ARIA® **POSTURAL BACKREST** is contoured to provide precise orientation within the wheelchair for optimal postural support. The foam is oversized to provide extra comfort and protection. The ARIA® **POSTURAL BACKREST** includes a movable foam lumbar pad that can be installed/inserted behind the existing foam cushion to provide additional positioning capability.

OUTER COVER

The outer cover is made of a breathable spacer mesh, combined with a leatherette material for added durability. Regular cleaning and inspection of the outer cover is recommended. Refer to Section 5.0 -Care and Maintenance.

2. SAFETY

IMPORTANT!

Check ALL parts for shipping damage. If shipping damage is noted, DO NOT use. Contact carrier/dealer for further instruction.

2.1 CUSTOMER TRANSPORT ON A WHEELCHAIR EQUIPPED WITH AN ARIA® BACKREST IN A MOTOR VEHICLE

The following guidelines should be made available to all parties who are responsible for the transportation of the client, such as schools and transport providers.

Wherever practicable Aria Wheels Srls recommends that the client is transferred out of the wheelchair into a vehicle seat using a lap and diagonal seat belt. ARIA® **POSTURAL BACKRESTS** may only be used in conjunction with wheelchairs that comply with the performance requirements of RESNA WC19 or ISO7176-19. For safe transportation of a client within a wheelchair in a motorvehicle, refer to the user manuals of the wheelchair and tie-down and occupant restraint system (WTORS).

Wheelchair Tie-down - The wheelchair must be securely tied down within the vehicle, in a forward facing position, in accordance with the wheelchair and tiedown manufacturers instructions.

Occupant Restraint - A suitable occupant restraint system must also be fitted in accordance with the manufacturer's instructions. A minimum is a lap and diagonal belt (lap belts on their own are not suitable).

Restraints that loop over the shoulder and anchor to the vehicle floor should be avoided where possible, as they can cause a heavy downward load through the client during an impact, the preferred type anchor above and behind the shoulder as with a car seat belt. The lap section of the belt should fit snugly over the pelvis of the client and should not be able to ride up to the abdominal area. The upper torso section of belt should be in contact with the client's chest and fit over the shoulder(s), while not cutting into the neck or slipping off the shoulder(s).

IMPORTANT !

If the wheelchair and ARIA® **BACKRESTS** are to be transported unoccupied, the wheelchair should be tied down with suitable equipment according to the manufacturer's instructions. If the equipment is to be dismantled for transportation, ensure that all parts are secured safely within the vehicle.

2.2 INTENDED USE

Our ARIA® **POSTURAL BACKRESTS** are designed to provide precise orientation within the wheelchair for optimal postural support

ATTENTION !

Do not use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as owner's manuals, service manuals or instruction sheets supplied with this product or optional equipment. If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, dealer or technical personnel before attempting to use this equipment - otherwise, injury or damage may occur.

2.3 ACCESSORIES WARNING

ARIA® products are specifically designed and manufactured for use in conjunction with ARIA® accessories. Accessories designed by other manufacturers have not been tested by Aria Wheels Srls and are not recommended for use with ARIA® products.

2.4 INSTALLATION WARNINGS

The procedures in this manual should be performed by a qualified technician. Attaching hardware that is loosely secured could cause loss of stability resulting in serious injury or damage. After any adjustments, repair or service and before use, make sure that all attaching component parts are secure.

DO NOT install the ARIA® POSTURAL BACKREST assembly onto back canes with an outside diameter greater than 1" (25mm) or less than 3/4" (19mm). Otherwise, injury or damage may occur. The ARIA® POSTURAL BACKREST is not intended for power seating applications. The mounting position of the ARIA® POSTURAL BACKREST is directly related to the wheelchair's stability. It may be necessary to reposition the front casters, rear wheels, back angle, tilt-in-space, recline position and/or seat depth before use. Ensure the ARIA® POSTURAL BACKREST is properly secured to the wheelchair before using. Otherwise injury or damage may occur.

ATTENTION !

The best way to avoid problems related to pressure sores is to understand their causes and your role in a skin management program. Your therapist and physician should be consulted if you have questions regarding individual limitations and needs. All back cushions should be selected carefully. Working with your therapist and physician is the best way to assure that a cushion choice matches your individual needs. As the needs of the individual become more complex, the evaluation becomes more important.

3. CONFIGURATION


3.1 HARDWARE KIT CONTENTS

IMPORTANT !

The screws are sold separately.
If any of the hardware is missing or misplaced, please contact our customer service department and arrangements will be made to send you the necessary replacements.


1. Fireproof foam padding with breathable coating (size chosen when ordering).
2. Rubber gasket for carbon fiber elements cod.VAR2
3. Stainless steel screws with cylindrical head, ultra-lowered, hexalobular hollow M6 x 16mm cod.BULL18
4. Stainless steel washer, hole diameter 6mm, external diameter 12mm cod.BULL22
5. Shell in carbon fiber or aluminum (size chosen when ordering)
6. Stainless steel screw with cylindrical head hexagon socket M6x16mm cod.BULL31
7. Universal clamp for wheelchair with rigid frame cod.ALCMT36 (composed of: ALCMT36.3_ Hardened pin d.6x50mm DIN 6325; + possible POM-C adapter for tube diameter 19mm ALCMT36.1 or ALCMT36.2 for tube diameter 22mm)

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3.3 MOUNTING CONFIGURATION

The procedures indicated in this manual must be performed by qualified technicians. Ensure the Matrix PB is properly secured to the wheelchair before using. Otherwise injury or damage may occur.

ATTENTION !

The installation hardware provided is high strength and tested for durability. **DO NOT** substitute hardware. Use only the hardware supplied.

1. Remove existing wheelchair back upholstery and tension back.
2. Install the universal clamp (ALCMT36) by lightly tightening the screws (BULL31) on the transverse tube of the wheelchair.
3. Place the shell as desired height by matching the holes on the back with the holes on the clamp.
4. Place the washers (BULL22) and tighten the supplied hexagon socket screws (BULL18).
5. Adjust the angle of the backrest as desired by tightening the hexagonal screws of the clamp (BULL31).
6. Position the padding paying attention to cover the edges of the backrest shell to avoid injury to the user.
7. In order to better determine the correct depth and the angle of the backrest, the user must be seated on the wheelchair.
8. The clamps and screws of the kit must be fully tightened after installing the backrest and after each adjustment and maintenance.

WARNING !

Improper configuration of the fixing pins can damage the back shell and/or the fixing clamp. Make sure that the pins are fully tightened; otherwise you could damage the back shell and/or the fixing clamp. **NEVER** use glue on the fixing pins. It is preferable after any repairs or adjustments to use the thread locker to prevent the pins from loosening over time.

4. USAGE

4.1 LUMBAR SUPPORT INSTALLATION/ADJUSTMENT

IMPORTANT !

The Lumbar Support is a pre-fabricated foam insert that provides lower back support for additional comfort. The position of the lumbar insert can be adjusted inside the ARIA® POSTURAL BACKREST cover, or it may be removed if no lumbar support is desired.

4.1 LUMBAR SUPPORT ADJUSTMENT

1. To access the lumbar support, open the hook and loop fastening strips on the shell.
2. Open the central zipper positioned on the back of the padding.
3. Put the foam lumbar support inside the Aria backrest cover between the foam back cushion and the back shell.
4. The lumbar support can be adjusted to any desired height/position.
5. Close the zipper on the back of the padding and reposition on the shell, making sure that it covers the edges to avoid injury to the user.

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5.2 INSPECTION

ATTENTION !

Inspect ALL fasteners weekly to ensure that mechanical connections and attaching hardware are tightened securely. DO NOT continue to use this product if problems are discovered. Corrective measures can be performed at, or arranged through your equipment supplier.

Perform a visual inspection of parts including hardware, brackets, upholstery materials, foams, and plastics for deformation, corrosion, breakage, wear or compression.

6. AFTER USE

6.1 REUSE

ATTENTION !

Danger of injury

Always have a dealer check that the product is not damaged before delivering the product to another user. If damage has been found DO NOT use the product. Otherwise, there is a risk of injury.

The back cushion is suitable for re-use.

The following actions must be taken:

- Cleaning and disinfection in accordance with the Cleaning 5.1.
- Inspection in accordance with Inspection Section 5.2.
- The back cushion should be adapted to the new user following the indications of the section Mounting the backrest 3.3 and Adjusting the lumbar support 4.1

6.2 DISPOSAL

Defend the environment providing an adequate disposal of your wheelchair. Observe the national law and local regulations for the disposal.

For a correct disposal, contact your specialized dealer or your council's public administration as to get the address of a local disposal centre.

7. TECHNICAL DATA

7.1 WEIGHT AND DIMENSIONS

All weight and dimension specifications refer to a width and height in all the possible configurations of the backrest. dimensions and weight may change according to the different configurations.

<i>Backrest width:</i>	<i>335 or 395 mm</i>
<i>Backrest height:</i>	<i>220 / 280 / 340 mm (height adjustable)</i>
<i>Backrest angle:</i>	<i>85° to 95° - adjustable respect to the floor</i>

<i>Backrests weight:</i>	<i>from 650 to 1000 gr*</i>
<i>Nominal load:</i>	<i>120 kg</i>



* The weight depends on the possible configurations

7.2 ENVIRONMENTAL CONDITIONS

Do not expose the backrest to temperatures below -20 °C and above 40 °C.

7.3 LABELS

The symbols mentioned on the label are those compliant with UNI EN ISO 15223-1:2012 and UNI CEI EN 1041:2013 regulations.

SIMBOLS	MEANING
	Catalogue number
	Manufacturer
	Date of Manufacture
	Read User Manual
SN	Serial number
CE	Compliant with CE regulations
	Rated flow
	Warning

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